



A β ₁₋₄₂, t-tau and p-tau₁₈₁ measurements in 120 PPMI CSF samples using xMAP/Luminex multiplex immunoassay

Ju Hee Kang, Teresa Waligorska, John Q Trojanowski and Leslie M Shaw, Department of Pathology & Laboratory Medicine and Institute on Aging, Center for Neurodegenerative Disease Research, Perelman School of Medicine, University of Pennsylvania

Summary

One hundred and twenty CSF samples collected according to the Parkinson's Progression Marker Initiative (PPMI) Research Biomarkers Laboratory Manual were analyzed for the measurement of three CSF biomarkers, amyloid beta (1-42) (A β ₁₋₄₂), total tau (t-tau) and tau protein phosphorylated at the threonine 181 position (p-tau₁₈₁) using the research-use-only multiplex xMAP Luminex platform and Innogenetics immunoassay kits (Innogenetics/Fujirebio, Ghent, Belgium). The never before thawed frozen aliquots of CSF samples stored at the PPMI Biorepository Core laboratory (Coriell Institute) were transferred to the University of Pennsylvania PPMI Bioanalytics core and stored at -80 °C until the day of analysis. Each study subject CSF sample aliquot was labeled with a unique code number. Using aqueous buffer-based standards of three biomarkers, calibration curves were generated for each biomarker. All standards, aqueous controls and CSF samples (including 2 CSF pools for assessment of run quality control) were analyzed in duplicate in each run (total of four runs for completion of the analyses). The % CV for replicates for the three biomarkers in two aqueous controls and two CSF pools are summarized in the table below. The %CV values for the duplicate analyses for PPMI study subject CSF samples were below 25%, except for t-tau (6 samples) showed over 25 % CV (26.9 – 29.0 %), and A β ₁₋₄₂ and p-tau₁₈₁ in 1 sample and t-tau in 2 samples showed low bead counts (< 50). These data are included in the report but we point these out as the results in these cases are not as precise as the remainder of the data.

Method

For each study subject, CSF was collected in a polypropylene syringe. The first 1-2 mLs of CSF was sent to a local laboratory for routine diagnostic testing and analyzed within 4 hours of collection for cell count, total protein and glucose as described in the Laboratory Manual. The remaining 15-20 mL of CSF was transferred into polypropylene tubes at room temperature, mixed gently, centrifuged at 2000×g for 10 minutes at room temperature and aliquoted into pre-cooled labeled polypropylene sample storage tubes followed by immediate freezing on dry ice. The frozen aliquots of CSF, n=120, were shipped overnight from Coriell Institute for Medical Research on dry ice to the PPMI Bioanalytics Core laboratory at the Perelman School of Medicine, University of Pennsylvania. Each aliquot was identified with a unique code number. A β ₁₋₄₂, t-tau, and p-tau₁₈₁ were measured in each of the never before thawed 120 CSF PPMI



aliquots using the research-use-only multiplex xMAP Luminex platform (Luminex Corp, Austin, TX) with Innogenetics INNO-BIA AlzBio3 (Ghent, Belgium) immunoassay kit-based reagents, as described previously (Olsson A et al, 2005; Shaw LM et al., 2011). In brief, kit reagents included well-characterized capture monoclonal antibodies specific for A β_{1-42} (4D7A3), t-tau (AT120), and p-tau₁₈₁ (AT270) chemically bonded to unique sets of color-coded microbeads, and analyte-specific detector antibodies (HT7 for t-tau and p-tau₁₈₁ and 3D6 for A β_{1-42}). Using aqueous buffer-based standards of three biomarkers at concentrations ranging from 22 to 2231 pg/mL for synthetic A β_{1-42} , 9 to 1359 pg/mL for recombinant tau, and 4 to 231 pg/mL for a tau synthetic peptide phosphorylated at the threonine 181 position (p-tau₁₈₁), calibration curves were generated for each biomarker. All standards, controls and CSF samples (including 2 CSF pools for run validation) were analyzed in duplicate. A 5 parameter logistic model with weighting(1/Y²) was used for calibration curve fitting. The lot numbers for the AlzBio3 immunoassay kits and calibration standards used in this study were: 215369 and 213775, respectively.

Table 1. QC data for the PPMI CSF biomarkers study

Aqueous Control		Run1	Run2	Run3	Run4	Mean \pm S.D.	Mean %CV
A β_{1-42}	A	159	170	161	157	162 \pm 5.8	3.60
	B	75	74	69	76	73 \pm 2.9	3.99
t-tau	A	103	90	94	87	93 \pm 7.0	7.45
	B	192	202	197	185	194 \pm 7.3	3.77
p-tau ₁₈₁	A	34	35	33	33	34 \pm 1.2	3.67
	B	44	47	42	44	44 \pm 2.3	5.17
CSF pools		Run1	Run2	Run3	Run4	Mean \pm S.D.	Mean %CV
A β_{1-42}	C01	187	216	215	179	195 \pm 18.9	9.67
	D11	199	181	205	189	193 \pm 10.4	5.37
t-tau	C01	153	170	158	158	159 \pm 8.1	5.09
	D11	52	50	54	45	51 \pm 3.9	7.69
p-tau ₁₈₁	C01	23	24	23	23	23 \pm 0.8	3.29
	D11	15	16	16	16	16 \pm 0.4	2.72



References

1. Olsson A, Vanderstichele H, Andreasen N, et al. Simultaneous measurement of β -amyloid(1-42), total tau and phosphorylated tau (thr181) in cerebrospinal fluid by the xMAP technology. Clin Chem 2005;51:336 –345
2. Shaw LM, Vanderstichele H, Knapik-Czajka M, et al. Qualification of the analytical and clinical performance of CSF biomarker analyses in ADNI. Acta Neuropathol 2011;121:597-609.
3. Parkinson's Progression Marker Initiative Research Biomarkers Laboratory Manual (Biologic Manual). <http://www.ppmi-info.org/wp-content/uploads/2011/05/PPMI-Biologics-Manual-April-2011-FINAL.pdf>

About the Authors

This document was prepared by Ju Hee Kang, Teresa Waligorska, John Q Trojanowski and Leslie M Shaw, Department of Pathology & Laboratory Medicine, Perelman School of Medicine, University of Pennsylvania. For more information please contact Ju Hee Kang at 215-662-6266 or by email at Ju.Kang@uphs.upenn.edu.

Notice: This document is presented by the author(s) as a service to PPMI data users. However, users should be aware that no formal review process has vetted this document and that PPMI cannot guarantee the accuracy or utility of this document.